



504

5 Section Couch

Providing additional comfort for prone lying patients, this 5 section is based on our popular Model 503 with the addition of a breathing hole and adjustable arm rests. Ideally suited to manipulation and examination and can have a wide range of accessories fitted.

This couch is also available with drainage facility as model 505.

Features: height adjustable arm rests, unique lifetime guarantee, easy operation, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

STANDARD SPECIFICATION


- Either electric or hydraulic height adjustment available.
- Able to lift up to 225kg (35 stone) evenly distributed from its lowest height.
- Electric models operated as standard by a hand switch (optional footswitch).
- Foot section adjustable between -20 and +90°.
- Backrest adjustment from -20 to +35°.
- Hoist clearance of 13cm (5").
- Adjustable foot for uneven floors.
- Fully washable/ fire retardant Ambla vinyl in a range of colours.
- Robust retractable castors for excellent mobility and safety.
- Couch goes horizontal for patient recovery in emergency situations.
- CE marked in conformity with MDD 93/42/EEC Annex IX Class 1.
- Electrical standard IEC601 part 1.

Standard Options

- Footswitch
- Battery back up (electric versions only)
- Paper roll holder (foot end only)

504 DIMENSIONS

Length: 185 cm (73")
Width: 63 cm (25")
Height: 46 cm (18") approx
to 95 cm (37") approx



504 Electric
504 Hydraulic

PACKED WEIGHTS AND DIMENSIONS FOR SHIPPING

Packed Weight: 85 Kgs
Shipping Dimensions: 190 x 66 x 64 cm

DELIVERY PROMISE

Our normal delivery service is 7 – 10 days from receipt of order; however if you have a specific requirement we will always endeavour to meet it.

FOR YOUR PEACE OF MIND

Our unique guarantee is for a LIFETIME on the mechanism and UNCONDITIONAL for 3 years on all components.. This even applies to accidental damage. We will not quibble. The only exclusion is damage to upholstery.

All of our products are CE marked. We are an ISO9001 audited company. Should you require any help concerning the most suitable product for your needs, we would be delighted to discuss it with you.

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000



VINYL

Upholstery

We use only the best quality vinyl proven over many years and ALL feature an anti microbial additive to help prevent cross infection.

All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. ALL include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products.

Note: Unless specified BEIGE will be supplied.



Beige



Grey



Black



Fawn



Chocolate



Brow



Damson



Wine



Ruby



Terracotta



Night




Ocean



Sky



Slate




Green



Sage



Lime



White

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

