

Overview:

Based on the popular model 503, this tilting version offers all the same features with the additional benefit of a trendelenburg seat tilt which offers greater comfort for patients and better positioning for the practitioner. Ideally suited for minor procedures, general examination and chiropody use.

Standard Features


- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Electric operation with trendelenburg tilt
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Hoist clearance 130mm



Overall Dimensions

Length: 186cm (73") approx
 Width: 63cm (25") approx
 Height: 48cm - 97cm approx
 (19" - 38") approx

Upholstery Dimensions



Recommended Options:



Double Foam




Foot Switch



D Cushion



Motor Back




Breathing Hole



Paper Roll Holder

Operational Parameters



Model 503T - Electric version

For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.

