

SEIRIN

Acupuncture Needles


G-type Needle



G-type


LONG NEEDLE WITH GUIDE TUBE

Using only the most refined materials and state-of-the-art technology, We created a premium acupuncture needle for the most effective treatment.




Needle tip is the newly developed SP (Seirin Point)

The latest grinding technology enabled the Seirin Point (SP), featuring a rounded tip.



The longer 25mm handle

The handle length is suitably proportioned for this long needle.



The package is designed for easy opening

The new blister package features "7mm seal top", "curved line" and "Embossed sheet".

PRODUCT LINEUP



No	5	8
C.CODE	Violet	Brown
Dia(mm)	0.25	0.30
Length75mm	○	○
Length90mm	○	○

